V01: Analize okoljske DNA


(angl. environmental DNA)

Laboratorijski praktikum
2026/2027

Asist. Aja Bončina


OSNOVE OKOLJSKE DNA (eDNA)


Starodavna DNA
(aDNA; ancient DNA)

Okoljska DNA *(eDNA; environmental DNA)*

I. GENETSKI MONITORING

nevretenčarska DNA *(iDNA; invertebrate DNA)*

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DNA Extraction and Purification from 250–1,000 mg Soil

LYSE

INHIBITOR REMOVAL

1. Transfer the glass beads from a Disruptor Tube to a 15 mL centrifuge tube. Add 0.25-1 g soil sample and 1 mL SLX-Mlus Buffer. Vortex at maximum speed for 3-5 minutes to lyse samples.

2. Add 100 µL DS Buffer. Vortex to mix thoroughly. Incubate at 70°C for 10 minutes. Briefly vortex the tube once during incubation. Centrifuge at 3,000 rpm for 3 minutes at room temperature.

3. Transfer 800 µL supernatant into a new 2 mL microcentrifuge tube (not provided). Add 270 µL chilled P2 Buffer. Vortex to mix thoroughly. Let sit on ice for 5 minutes. Centrifuge at maximum speed for 5 minutes.

4. Carefully transfer the supernatant to a new 2 mL microcentrifuge tube. Add 0.7 volumes 100% isopropanol. Mix thoroughly by inverting tube for 20-30 times. If the soil contains very low DNA, incubate the sample at -20°C for 1 hour. Centrifuge at maximum speed for 10 minutes. Carefully aspirate and discard the supernatant. Do not disturb the DNA pellet.

5. Invert the tube on absorbent paper for 1 minute to drain the liquid. Add 200 µL Elution Buffer. Vortex for 10 seconds. Incubate at 70°C for 10-20 minutes to dissolve the DNA pellet.

6. Add 100 µL cHTR Reagent that has been completely resuspended. Vortex to mix thoroughly. Let sit at room temperature for 2 minutes. Centrifuge at maximum speed for 2 minutes.

7. Transfer the cleared supernatant to a new 2 mL microcentrifuge tube. If supernatant still has a dark color from the soil, repeat Step 6 for a second cHTR Reagent step. This will require additional cHTR Reagent (Cat# CHTR-50) to be purchased separately.

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BIND

8. Add an equal volume of XP1 Buffer. Vortex to mix thoroughly.

9. Insert a HiBind® DNA Mini Column into a 2 mL Collection Tube. Transfer the sample from Step 8 to the HiBind® DNA Mini Column. Centrifuge at 10,000*g* for 1 minute at room temperature. Discard the filtrate and reuse the Collection Tube.

WASH

10. Add 500 µL HBC Buffer diluted with 100% isopropanol. Centrifuge at 10,000*g* for 1 minute. Discard the filtrate and the Collection Tube.

11. Transfer the HiBind® DNA Mini Column into a new 2 mL Collection Tube. Add 700 µL DNA Wash Buffer diluted with 100% ethanol. Centrifuge at 10,000*g* for 1 minute. Discard the filtrate and reuse the Collection Tube.

12. Centrifuge the empty HiBind® DNA Mini Column at maximum speed for 2 minutes at room temperature. This step is critical in removing residual ethanol that may interfere with downstream applications.

ELUTE

13. Transfer the HiBind® DNA Mini Column into a clean 1.5 mL microcentrifuge tube. Add 50-100 µL Elution Buffer heated to 70°C directly onto the center of HiBind® matrix. Let sit at room temperature for 1-2 minutes. Centrifuge at maximum speed for 1 minute.

14. Take the filtrate from Step 13 and place onto the center of the same HiBind® DNA Mini Column used in the procedure. Let sit at room temperature for 1 minute. Centrifuge at maximum speed for 1 minute.

15. Store eluted DNA at -20°C.